



- ◆ High performance
- ◆ Widest gradient range
- ◆ Single-screen programming
- ◆ Flexible platform and consumables
- ◆ Space-saving design

### Amplitronyx Thermal Cyclers

Bring cutting-edge design and exceptional flexibility to your laboratory with high performance Amplitronyx Thermal Cyclers!

Compact and powerful, Amplitronyx cyclers are incredibly easy to use. Large touch-screens allow you to input your program on a single screen, and thermal blocks can be easily interchanged without tools.

| ←     | Run   | Program | Utilities     | Block        |
|-------|-------|---------|---------------|--------------|
| Stage | Cycle | Step    | Temp. Options | Time Options |
| 1     | 1     | 1       | 94.0          | 1:00         |
| 2     | 30    | 1       | 94.0          | 0:30         |
|       |       | 2       | 50-70         | 0:30         |
|       |       | 3       | 72.0          | 0:30         |
| 3     | 1     | 1       | 72.0          | 5:00         |
|       |       | 2       | 4.0           | HOLD         |

Program Pref. Add Delete Print Save as  
User I Standard I Editing

Amplification reactions take place under tightly regulated thermal conditions where samples are guaranteed to stay at the programmed temperature for the entire length of the step. The heated lid minimizes condensation within the sample tube, and the thermal block interior is protected from condensation by unique double-seal construction.


Amplitronyx thermal blocks are engineered with precise attention to well size, shape and placement. Meticulous design and lightweight construction results in faster heating and cooling than that found in standard blocks. Many years of engineering and application expertise ensures extraordinary block uniformity and reproducible performance. Cyclers are backed by first-class technical support and a two-year warranty.

#### Amplitronyx 4

The Amplitronyx 4 is the most cost-effective option for a powerful, high-performance cycler. Store up to 100 seven-stage/step programs for touchdown and standard amplification in either 0.2 ml or 0.5 ml tubes.

#### Amplitronyx 6

If experiments require a temperature gradient or the use of multi-well plates, the Amplitronyx 6 offers the widest gradient range available and allows for tool-free swapping of four different thermal blocks. Program a 36°C-wide gradient across the 0.2 ml and 0.5 ml blocks, or increase your throughput and experimental flexibility even more with 384-well and slide blocks. The temperature of the heated lid can be adjusted to suit your protocol, and up to 200 ten-stage/step programs for touchdown and standard amplification can be stored.



## Specifications

| Sample Blocks                          | Amplitrionyx™ 4   | Amplitrionyx™ 6   |
|--|---|---|
| Capacity                               | <b>CombiBlock:</b> 80 x 0.2 ml tubes or 32 x 0.5 ml tubes | <b>0.2 ml Gradient Block:</b> 96 x 0.2 ml tubes or 1x 96-well plate<br><b>0.5 ml Gradient Block:</b> 54 x 0.5 ml tubes<br><b>384-Well Block:</b> 1x 384-well plate<br><b>Flat Block:</b> 4 slides |
| Temperature                            | 4°C – 99°C  |   |
| Heating rate                           | Up to 3°C /sec  |   |
| Cooling rate                           | Up to 2°C /sec  |   |
| Uniformity                             | ±0.5°C within 15 sec                                      | ±0.4°C within 15 sec  |
| Material                               | Nickel-plated aluminum                                    |   |
| <b>Heated Lid</b>                      |   |   |
| Temperature                            | 105°C (fixed)   | 99 – 105°C (adjustable)   |
| Pressure                               | Spring-loaded   | Electromechanical   |
| <b>Programming</b>                     |   |   |
| Display/Input                          | 5.7" diag. touchscreen (b&w)                              | 5.7" diag. touchscreen (color)  |
| Number of programs                     | 100   | 200   |
| Maximum program stages                 | 7   | 10  |
| Maximum steps per stage                | 7   | 10  |
| Gradient software                      | No  | 36° range from 30°C– 80°C   |
| Max programmed dwell time              | Infinite  |   |
| Pause function                         | Yes   |   |
| Temp increments/decrements (Touchdown) | Yes   |   |
| Time increments                        | Yes   |   |
| Auto-restart function                  | Yes   |   |
| Run "end time" calculation             | Yes   |   |
| Temperature control options            | Block control and simulated tube control                  | Block control, active tube control, and simulated tube, plate and slide controls  |
| <b>Working Conditions</b>              |   |   |
| Temperature                            | Ambient (12°C – 32°C)                                     |   |
| Relative humidity                      | Up to 80%   |   |
| Power requirements                     | 533W at 115/230VAC +/- 10% and 50/60Hz                    |   |
| <b>Physical Characteristics</b>        |   |   |
| Cycler Dimensions (W x D x H)          | 10 x 10 3/4 x 10 inches (250 x 274 x 250 mm)              |   |
| Shipping Dimensions                    | 13 1/2 x 13 1/2 x 12 inches (34 x 34 x 30 cm)             |   |
| Shipping Weight                        | 26 lbs (12 kg)  |   |

## Ordering Information

| Catalog # | Description   |
|-----------|---|
| A4C       | Amplitrionyx™ 4 Thermal Cycler with CombiBlock            |
| A62       | Amplitrionyx™ 6 Thermal Cycler with 0.2 ml Gradient Block |
| A65       | Amplitrionyx™ 6 Thermal Cycler with 0.5 ml Gradient Block |
| A63       | Amplitrionyx™ 6 Thermal Cycler with 384-well block        |
| A6F       | Amplitrionyx™ 6 Thermal Cycler with Flat Block            |
| A6-P      | Optional Printer for Amplitrionyx™ 6                      |
| A6-B2     | 0.2 ml Thermal Block for Amplitrionyx™ 6                  |
| A6-B5     | 0.5 ml Thermal Block for Amplitrionyx™ 6                  |
| A6-B3     | 384-Well Thermal Block for Amplitrionyx™ 6                |
| A6-BF     | Flat Block for Amplitrionyx™ 6                            |

Rev: 020108

